

# PHILIPS

Sonic electric  
toothbrush with app

DiamondClean  
9000

Built-in pressure sensor

Smart brush head recognition

4 modes, 3 intensities

Connected brushing made easy



HX9911/13



## Whiter, healthier teeth for life

Whiter teeth in just 1 day\*

Philips Sonicare's best whitening in our most elegant sonic electric toothbrush. Switch to Philips Sonicare.

### **DiamondClean users experience improved oral health**

- Remove up to 100% more stains\* in just 3 days
- Stains disappear for a brighter, whiter smile

### **Stay on track with a personalized progress report**

- Start and keep up healthy habits

### **Let the DiamondClean 9000 adapt to your needs**

- Personalized brushing experience
- Lets you know if you're pressing too hard
- BrushSync automatically selects the best mode for you
- Always know when to replace your brush heads

### **Enjoy a stylish and user-friendly design**

- Simple to charge and style to boot
- Charge for up to two weeks of battery life


# Highlights

## Remove up to 5x more stains\*\*\*




Get rid of surface stains and enjoy a whiter smile with the W3 Premium White brush head. The brush head's densely packed bristles are especially designed to tackle stains and it's clinically proven to remove up to 100% more stains in just three days.\*\*\*

## W3 Premium White brush head



Click on the W3 Premium White brush head to remove surface stains and reveal your whitest smile. Its densely packed central stain removal bristles whiten your teeth from day one.

## Progress Report



The DiamondClean 9000 gives you the guidance you need to improve and maintain healthy brushing habits in between dentist checkups. Built-in smart sensors let you know when you're using too much pressure, and by connecting your brushing experience to the Sonicare app, a personalized Progress Report helps you stay on track, to see how much you've improved over time.

## 4 modes, 3 intensities




The DiamondClean 9000 comes with Clean, White+, Gum Health and Deep Clean+ to take care of your brushing needs. Clean mode is for exceptional everyday cleaning, White+ is for removing stains, Gum Health provides a gentle yet effective clean for gums and Deep Clean+ gives you an invigorating deep clean. Three intensities allow you to choose between a higher setting to boost your clean and a lower one for more-sensitive mouths.

## BrushSync mode pairing



Smart brush heads ensure you're using the right mode and intensity for the best possible clean. For an example, say you're using the W3 Premium White brush head. Your DiamondClean 9000's BrushSync technology will automatically sync your brush head with the White+ mode to help whiten your teeth.

## Built-in pressure sensor



You might not notice if you're brushing too hard, but your DiamondClean 9000 will. If you need to ease up, the toothbrush will make a pulsing sound. It's a heads up to let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

## BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. The BrushSync replacement reminder on your handle and a short beep gives you a heads up when it's time to replace it.

## Long-life battery



Your DiamondClean 9000 will last for up to two weeks on a single charge with regular use.

## Design with storage in mind



Enjoy two weeks of regular use from a single full charge. Your deluxe USB travel case doubles as a charger so you stay topped up on the go. Place your toothbrush in the case and connect to your laptop or plug into a wall socket. Also includes a brush head holder for extra hygiene when traveling. As for charging at home, our sleek charger glass fits stylishly into your bathroom and doubles as a glass for rinsing after brushing. Simply place your brush into the glass.

# Specifications

## Modes

**3 intensities:** Low, Medium, High

**Clean:** For exceptional everyday clean

**Deep Clean+:** For an invigorating deep clean

**Gum Health:** Special attention to molars

**White+:** To remove surface stains

## Items included

**Handle:** 1 DiamondClean 9000

**Brush heads:** 1 W3 Premium White

**Travel case:** USB charging travel case

**Glass charger:** 1

**Charger base:** 1

## Connectivity

**Bluetooth® wireless technology:** Connected brushing app

## Compatibility

**Android compatibility:** Android phones, Bluetooth 4.0 enabled tablets

**iOS compatibility:** iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher, with iOS7 operational system

## Design and finishing

**Color:** Black

## Cleaning performance

**Whitening benefits:** Whitens teeth in 1 day\*

**Pressure feedback:** Vibrates handle to alert user

**Timer:** BrushPacer and SmarTimer

## Smart sensor technology

**Pressure sensor:** Alerts when brushing too hard

**BrushSync Replacement Reminder:** Always

know when to, replace brush heads

**BrushSync Technology:** smart brush head, Connects smart handle and

## Ease of use

**Handle compatibility:** Easy click-on brush heads

**Replacement reminder:** To always ensure best results, reminder icon lights up

**Handle:** Slim ergonomic design

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Operating time (full to empty):** 14 days\*\*\*

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

