

**NEW**

# TOP TECHNOLOGY MADE BY SCHICK



User-friendly & practice-oriented

### Outstanding performance

The new Q3 series from Schick offers unbeatable performance parameters with up to 60,000 rpm and 11 Ncm. Thanks to new technology, the speed can be controlled extremely sensitively and is therefore ideally suited for processing even the most demanding dental materials.

### Multifunctional and reliable

All important parameters are clearly displayed on the bright OLED display and are easy to read from all angles. Additional functions such as clockwise/anticlockwise rotation and cruise control leave nothing to be desired and can be conveniently activated via direct selection keys.

### Three power classes

With three different performance classes, the Q3 covers the requirements of a wide range of applications. The three unit types are supplied with the following handpieces as standard:

- ▶ The **Q3 Basic** with Q Basic motor handpiece, which is perfectly suited for universal use in everyday laboratory work thanks to its robust design.
- ▶ The **Q3 Profi** with Q Profi motor handpiece, which not only lies comfortably in the hand due to its soft-grip surface, but also scores with optimal noise and heat insulation.
- ▶ The **Q3 Premium** with Q3 motor handpiece – the ultra-light and ergonomically perfectly balanced motor handpiece, which with its robust inner workings is also suitable for the heaviest grinding, cutting and milling work.

**With Q3 you get absolute top technology Made by Schick.**



#### Q3 Basic

Knee control complete Art. No. 10765  
Foot control complete Art. No. 10768  
Table-top control complete Art. No. 10773

#### Q3 Profi

Knee control complete Art. No. 10775  
Foot control complete Art. No. 10778  
Table-top control complete Art. No. 10783

#### Q3 Premium

Knee control complete Art. No. 10785  
Foot control complete Art. No. 10788  
Table-top control complete Art. No. 10793



**Useful footpedals as additional control options:**  
Footswitch (on/off) Art. No. 6370/2  
Dynamic foot control Art. No. 9440

### Technical data

	Q3 Basic	Q3 Profi	Q3 Premium
Number of revolutions	200 – 40,000 rpm	200 – 50,000 rpm	200 – 60,000 rpm
Torque	7 Ncm	9 Ncm	11 Ncm
Operating voltage	100 – 230 Volt	100 – 230 Volt	100 – 230 Volt
Operating frequency	50/60 Hz	50/60 Hz	50/60 Hz
Power	100 Watt	150 Watt	200 Watt
Overload protection	No waiting time when restarting	No waiting time when restarting	No waiting time when restarting
Warranty	1 year 5 years electric motor	1 year 5 years electric motor	1 year 5 years electric motor

Handpiece			
Diameter	Max. 30 mm	Max. 28 mm	Max. 29 mm
Length	158 mm	158 mm	144 mm
Weight without cable	240 g	240 g	195 g
Weight with cable	370 g	370 g	330 g
Cooling	Self sealing system without cooling	Self sealing system without cooling	Self sealing system without cooling
Drive system	Brushless DC-Motor	Brushless DC-Motor	Brushless DC-Motor
Concentricity	< 0.02 mm	< 0.02 mm	< 0.02 mm
Chucks	Ø 2.35 mm (standard) Ø 3.00 mm (optional)	Ø 2.35 mm (standard) Ø 3.00 mm (optional)	Ø 2.35 mm (standard) Ø 3.00 mm (optional)
Tool change	Quick release	Quick release	Quick release

Control units			
Dimensions	Knee device	Foot device	Table-top device
Width	128 mm	175 mm	105 mm
Height	208 mm	95 mm	200 mm
Depth	207 mm	260 mm	210 mm
Weight	1,360 g	2,300 g	1,880 g

**Scope of delivery:**  
Control unit, handpiece with cable,  
handpiece rack, mains cable

NEW

## Powerstick

Unique solution for more performance

### More Power

With the newly developed Powerstick from Schick, the two power classes Q3 Basic and Q3 Profi can be retrofitted with additional power if required.

### Your benefits

- ▶ Increase in speed by 10,000 rpm and 2 Ncm torque.
- ▶ Simple and fast use of the Powerstick – no conversions required.
- ▶ Subsequently a powerful plus in performance without the need to purchase an additional micromotor.
- ▶ Not device-specific, can be used flexibly at different laboratory workplaces.

Powerstick  
Art. No. 10745



# FOR THE HIGHEST DEMANDS

# MADE BY SCHICK

THE BEST FOR YOUR LABORATORY



NEW



[www.schick-dental.de](http://www.schick-dental.de)



Schick GmbH  
Lehenkreuzweg 12  
D-88433 Schemmerhofen  
Telefon +49 7356 9500-0  
Telefax +49 7356 9500-95  
info@schick-dental.de  
www.schick-dental.de



Your local dealer

Subject to technical modification without prior notice

21793\_9/21 e

 **SCHICK**  
create the future