


Dalbo®-Rotex according to Dr. Dalla Bona.





Root canal anchor combined with prosthetic retainer.



Characteristics


- The insertion direction can deviate up to 6°
- Highly corrosion-resistant material: made of pure titanium
- Slightly tapered
- Cement flow-off grooves
- Rounded post end
- The construction can be prepared in a single session
- For cutting the thread, only 2 – 3 turns are necessary
- The cement flow-off grooves enable tension-free cementing

Image	Item	Description	Cat. No.
Complete parts			
	Dalbo®-Rotex GT (size 1)	Delivery unit: pack of 5 pieces	051 864
	Dalbo®-Rotex GT (size 2)	Delivery unit: pack of 5 pieces	051 865
Components			
	Female part G	Mouth-resistant plastic, can be integrated directly into the denture plastic Delivery unit: pack of 5 pieces	051 868
	Dalbo®-B female part E	Precious metal alloy, can be integrated directly into the denture plastic	051 511
	Male part T (size 1)	Pure titanium (Grade 4)	051 869
	Male part T (size 2)	Delivery unit: pack of 5 pieces	051 870
Auxiliary instruments			
	Root canal drill (white) Anchor size 1		070 380
	Root canal drill (yellow) Anchor size 2		070 381
	Face cutting drill (white) Anchor size 1		070 499
	Face cutting drill (yellow) Anchor size 2		072 474
	Reamer (white) Anchor size 1		070 298
	Reamer (yellow) Anchor size 2		070 299

Image	Item	Description	Cat. No.
	Spanner key for spherical head		070 498
	Thomas spanner key	for trimmer	070 221
	Coupling piece	For extending instruments, suitable for angled handpiece and Thomas spanner key	072 414
	Instrument for mounting of elastomeric ring	For Dalbo®-B female part E (051 511)	070 205
	Activator	For Dalbo®-B female part E (051 511)	070 197
	Deactivator	For Dalbo®-B female part E (051 511)	070 199
Auxiliary parts			
	Spacer G	Spacers replace the anchoring female parts during plastic polymerisation in the laboratory.	070 440
	Transfer jig	For fabricating the master model	070 511
Set			
	Introduction set Dalbo®-Rotex (Brunner/Dalla Bona)	Consisting of: 2 Dalbo®-Rotex GT according to Dr. Dalla Bona 051 864 2 Dalbo®-Rotex GT according to Dr. Dalla Bona 051 865 2 Dalbo®-Rotex GT according to Prof. Brunner 051 866 2 Dalbo®-Rotex GT according to Prof. Brunner 051 867 1 Root canal drill 070 380 1 Root canal drill 070 381 1 Face cutting drill 070 499 1 Face-cutting drill 072 474 1 Face-cutting drill 070 295 1 Face-cutting drill 072 473 1 Reamer 070 298 1 Reamer 070 299 1 Thomas spanner key 070 221 1 Spanner key 070 498 1 Spanner key 070 500 8 Spacers 070 440 2 Transfer jigs 070 511 2 Transfer jigs 070 512	051 904
	Endobox Rotex/Rotex-RD/Dalbo®-Rotex	Empty, for root canal anchors Dalbo®-Rotex, without anchors and instruments	08000002

Dalbo®-Rotex according to Prof. Brunner.






Root canal anchor combined with prosthetic retainer.



Characteristics

- The insertion direction can deviate up to 18°
- Highly corrosion-resistant material: made of pure titanium
- Slightly tapered
- Cement flow-off grooves
- Rounded post end
- The construction can be prepared in a single session
- For cutting the thread, only 2 – 3 turns are necessary
- The cement flow-off grooves enable tension-free cementing

Image	Item	Description	Cat. No.
Complete parts			
	Dalbo®-Rotex GT (size 1)	Delivery unit: pack of 5 pieces	051 866
	Dalbo®-Rotex GT (size 2)	Delivery unit: pack of 5 pieces	051 867
Components			
	Female part G	Mouth-resistant plastic, can be integrated directly into the denture plastic Delivery unit: pack of 5 pieces	051 868
	Dalbo®-B female part E	Precious metal alloy, can be integrated directly into the denture plastic	051 511
	Male part T (size 1)	Pure titanium (Grade 4)	051 871
	Male part T (size 2)	Delivery unit: pack of 5 pieces	051 872
Auxiliary instruments			
	Root canal drill (white) Anchor size 1		070 380
	Root canal drill (yellow) Anchor size 2		070 381
	Face cutting drill (white) Anchor size 1		070 295
	Face cutting drill (yellow) Anchor size 2		072 473
	Reamer (white) Anchor size 1		070 298
	Reamer (yellow) Anchor size 2		070 299

Image	Item	Description	Cat. No.
	Spanner key for spherical head		070 500
	Thomas spanner key	for trimmer	070 221
	Coupling piece	For extending instruments, suitable for angled handpiece and Thomas spanner key	072 414
	Instrument for mounting of elastomeric ring	For Dalbo®-B female part E (051 511)	070 205
	Activator	For Dalbo®-B female part E (051 511)	070 197
	Deactivator	For Dalbo®-B female part E (051 511)	070 199
Auxiliary parts			
	Spacer G	Spacers replace the anchoring female parts during plastic polymerisation in the laboratory.	070 440
	Transfer jig	For fabricating the master model	070 512
Set			
	Introduction set Dalbo®-Rotex (Brunner/Dalla Bona)	Consisting of: 2 Dalbo®-Rotex GT according to Dr. Dalla Bona 051 864 2 Dalbo®-Rotex GT according to Dr. Dalla Bona 051 865 2 Dalbo®-Rotex GT according to Prof. Brunner 051 866 2 Dalbo®-Rotex GT according to Prof. Brunner 051 867 1 Root canal drill 070 380 1 Root canal drill 070 381 1 Face cutting drill 070 499 1 Face-cutting drill 072 474 1 Face-cutting drill 070 295 1 Face-cutting drill 072 473 1 Reamer 070 298 1 Reamer 070 299 1 Thomas spanner key 070 221 1 Spanner key 070 498 1 Spanner key 070 500 8 Spacers 070 440 2 Transfer jigs 070 511 2 Transfer jigs 070 512	051 904
	Endobox Rotex/Rotex-RD/Dalbo®-Rotex	Empty, for root canal anchors Dalbo®-Rotex, without anchors and instruments	08000002